



Vacotron 460

Suction module:
A valuable addition
to the Endomed 482
or Sonopuls 492

Vacotron 460

The Vacotron 460 is a module for the 4 series. You can apply electrotherapy by using vacuum electrodes with this suction module. This is accomplished fast and simple. Positioning electrodes could not be any easier.

The Vacotron 460 is a valuable addition to the electrotherapy devices of the 4 series.

- Can be connected to all the electrotherapy and combination units of the 4-series
- Quick application of vacuum electrodes on positions where fixation straps cannot be used
- Continuous and intermittent suction
- Accurate adjustment of the suction strength
- Unique cable connection for fast application of rubber electrodes



The Vacotron 460 can only be used in combination with the Sonopuls 492 or the Endomed 482.



Vacotron 460


1498950 Vacotron 460 module

Can only be used in combination with the Sonopuls 492 or the Endomed 482

Standard accessories Vacotron 460

- 1498756 4-Series Information Booklet
- 2x 3444503 Vacuum electrodes Ø 60 mm, set of 2
- 3444505 Sponges Ø 65 mm, set of 4
(for vacuum electrodes Ø 60 mm)
- 2x 3444507 Vacuum lead hose red
- 2x 3444508 Vacuum lead hose black


ORDERING DATA



Vacotron 460

TECHNICAL SPECIFICATIONS

Dimensions (WxDxH)	24 x 28.6 x 9.3 cm
Weight	2 kg
Vacuum	continuous and pulsed, 0-800 mbar, continuously adjustable
Pulsed Vacuum	pulse : pause = 0.5 : 0.5 seconds, 1 : 1 second



Vacotron 460

ACCESSORIES

Vacuum accessories

- 3444509 Vacuum electrodes Ø 30 mm, set of 2
- 3444503 Vacuum electrodes Ø 60 mm, set of 2
- 3444504 Vacuum electrodes Ø 90 mm, set of 2
- 3444516 Sponges Ø 30 mm, set of 4 (for vacuum electrodes Ø 30 mm)
- 3444505 Sponges Ø 65 mm, set of 4 (for vacuum electrodes Ø 60 mm)
- 3444506 Sponges Ø 95 mm, set of 4 (for vacuum electrodes Ø 90 mm)
- 3444507 Vacuum lead hose, red
- 3444508 Vacuum lead hose, black



3444509
Ø 30 mm



3444503
Ø 60 mm



3444504
Ø 90 mm



3444516
Ø 30 mm



3444505
Ø 65 mm




3444506
Ø 95 mm



3444507



3444508



RELATED PRODUCTS

Sonopuls 492

The top device of the 4 series is the Sonopuls 492. This complete unit ensures that you can apply ultrasound therapy, electrotherapy and combination therapy.

Endomed 482

The Endomed 482 is the electrotherapy device for the demanding therapist.

EN-Car

The EN-Car is your companion for any piece of Enraf-Nonius equipment. For laser, ultrasound, biofeedback, shockwave or electrotherapy, with the EN-Car you have it all close at hand.



Please also see the
separate product leaflets at
www.enraf-nonus.com





Vacotron 460

Suction module:
A valuable addition
to the Endomed 482
or Sonopuls 492