



## Manumed Exercise couch

For safe and comfortable  
neurological and  
exercise therapy



PARTNER FOR LIFE

# Manumed Exercise couch

Neurological exercise therapy with the Manumed Bobath. The large, extra wide treatment surface of 100 cm to 120 cm makes it easy to support and assist patients during therapeutic exercises. The Bobath and Vojta motion methods can be applied to perfection owing to the unmatched stability. A sturdy couch affording both patient and practitioner the most comfortable treatment positions.

Neurological exercise therapy with the Manumed Bobath. The large, extra wide treatment surface of 100 cm to 120 cm makes it easy to support and assist patients during therapeutic exercises. The Bobath and Vojta motion methods can be applied to perfection owing to the unmatched stability. A sturdy couch affording both patient and practitioner the most







## Manumed Exercise couch


- Very large surface (100 or 120x200 cm)
- Wheel raising mechanism for easy movement of the couch
- Electric or hydraulic height adjustment
- The 2-section model can be used as a neurological exercise or treatment couch




# Manumed Exercise couch 2-section height adjustable

Type

- 511 Manumed Exercise 2-section, hydraulic H/L
- 521 Manumed Exercise 2-section, electric H/L



## ORDERING DATA






# Manumed Exercise couch 1-section height adjustable

Type

- 510 Manumed Exercise 1-section, hydraulic H/L
- 520 Manumed Exercise 1-section, electric H/L



## ORDERING DATA



# Manumed Exercise couch


## TECHNICAL SPECIFICATIONS

Size tabletop (lxw):	200 x 100 cm for models 510 and 520 optional 120 cm width for models 510 and 520 203 x 120 cm for models 511 and 521
Height adjustment:	44-97 cm
Lifting time (min.-max.):	approximately 18 sec.
Lifting capacity:	250 kg for models 520 and 521 (electric) 200 kg for models 510 and 511 (hydraulic)
Power supply:	120/230VAC 50/60 Hz
Current consumption:	2.0 A max.
NOTE:	ALL-AROUND OPERATION IS AN OPTION.

This is a medical device



It complies with all applicable medical regulations. For more details please consult the instructions for use (downloadable via [www.enraf-nonus.com](http://www.enraf-nonus.com)).





# Manumed Colours

## OPTIONS<sup>(1)</sup>

A matching Enraf-Nonius couch for every interior. As in other practices, a pleasant ambience is important in physiotherapy practices. Clients must feel comfortable in the room where they are being treated. As each practice is furnished differently and has its own colour scheme, we regard it as important for the treatment couches to fit in with your look and feel. That is why we can manufacture our treatment couches in several frame and upholstery colours. Our imitation leather upholstery contains no phthalates and is flame retardant – a unique combination that makes this material suitable for medical use!

### NATURAL

- 811 volcano black
- 850 walnut
- 851 sand
- 886 leaf
- 888 snow


### LIFESTYLE

- 410 anthracit
- 406 petrol
- 407 cherry
- 408 terracotta
- 409 skin

### EXECUTIVE

- 012 black
- 213 grey
- 199 lavender
- 204 sirius blue
- 169 cornflower blue
- 170 baby blue


(1) = NOTE: Please order options together with the couch because of modification of the couch during production.





# Manumed Colours

## NATURAL



811 volcano black



850 walnut



851 sand



886 leaf



888 snow



# Manumed Colours


## LIFESTYLE



410 anthracit



406 petrol



407 cherry



408 terracotta




409 skin



# Manumed Colours


## EXECUTIVE




012 black




213 grey



199 lavender



204 sirius blue



169 cornflower blue



170 baby blue



## Frame colours

Ral 7035 grey

Ral 9010 white

OPTIONS<sup>(1)</sup>



NOTE: Please order options together with the couch because of modification of the couch during production.



# Manumed Exercise couch

## OPTIONS<sup>(1)</sup>


### Upholstery options

Section Width	100 cm	120 cm
Standard 4 cm thick	s	s
510 / 520	s	o
511 / 521	x	s


s = standard 0 = optional x = not available



100 cm



120 cm



Standard

### Foam density:

Standard: 38 Kg/m<sup>3</sup>

### Height adjustment

- A. Hydraulic (510/511)
- B. Electric (520/521) Footswitch (standard)
- C. Electric (520/521) All-around (optional)



A



B



C



(1) = NOTE: Please order options together with the couch because of modification of the couch during production



# Manumed Exercise couch

## OPTIONS<sup>(1)</sup>

### Lifting capacity

A. Electric - max. 250 kg (standard 520/521)

B. Hydraulic max 200 kg (510/511)



A



B

### Optional features

Safety switch



### Fixation straps


3445474 Fixation strap with safety lock 350 cm



3445474

## ACCESSORIES

(1) = NOTE: Please order options together with the couch because of modification of the couch during production



# Stools

- 3xxx010 Stool, height adjustable 40-53 cm
- 3xxx011 Stool, mobile, height adjustable 59-85 cm
- 3xxx012 Stool, mobile, provided with footrest, height adjustable 59-85 cm
- 3xxx013 Stool, mobile, anatomic seat, height adjustable 43-56 cm
- 3xxx020 Stool, mobile, saddle seat, height adjustable 59-85 cm



3xxx010



3xxx011



3xxx012













3xxx013




3xxx020

Available in all colours shown below

NATURAL	LIFESTYLE	EXECUTIVE
 <span>811</span>  <span>851</span>	 <span>410</span>  <span>408</span>	 <span>012</span>  <span>213</span>  <span>199</span>  <span>204</span>  <span>169</span>  <span>170</span>


# ACCESSORIES





## Rolls and cushions

3xxx002	Halfroll, flat bottom, small, especially for the ankle	28x9x11 cm
3xxx003	Halfroll, flat bottom, large	56x16x18 cm
3xxx023	Massage roll	50x Ø15 cm
3xxx004	Neck cushion, hourglass type	30x15 cm
3xxx005	Massage cushion	30x28x8 cm
3xxx015	Wedge shaped cushion	50x30x15 cm



Available in all colours shown below

NATURAL	LIFESTYLE	EXECUTIVE

## ACCESSORIES



## Rolls and cushions

1550813 Plug for nosehole, triangular shaped <sup>(1)</sup>


16,5 x 21 cm

3445297 Wedge shaped, Kaltenborn, rubber model <sup>(1)</sup>

22,5 x 10,5 x 8,5 cm

3446161 Support bag with sand 1,2 kg <sup>(2)</sup>

15x25 cm



1550813



3445297




3446161

(1) = only available in black

(2) = only available in grey

## ACCESSORIES







# Manumed Exercise couch

For safe and comfortable  
neurological and  
exercise therapy

Enraf-Nonius B.V. | Vareseweg 127 | 3047 AT Rotterdam | The Netherlands  
www.enraf-nonius.com | info@enraf-nonius.nl | + 31 - ( 0 ) 10 20 30 600



PARTNER FOR LIFE