


PHOTON LASER

CLASS IV




HEAVENLY LIGHT






Quick Treatment Selection



Inbuilt Pathology Library



Patient Attributes

Photon Laser Options

- 10 Watts - 980 nm (33000)
- 10 Watts - 1064 nm (33001)
- 10 Watts - 980/810 nm (33002)
- 15 Watts - 1064 nm. (33003)

Features

- Continuous
- Pulsed: Frequency Variable up to 2000Hz
- Duty Cycle variable up to 50%
- High-Intensity Pulse: Sharp pulses for treatment at maximum power and depth with adjustable heat feeling
- Inbuilt Pathology Library
- Large 7 inch Colour Touch Screen

Standard Accessories




Laser Applicator




Spacer






Wave Length Mixture



Wave Settings



Treatment Ready Screen

Standard Accessories With Trishul



Laser Applicator



Spacer



Joint/ Point
Applicator



Spherical/ Massage
Applicator

And CART

1064 nm

810 nm

980 nm



TRISHUL

A Laser with Three Wavelength

Three different interactivity
within biological tissue

Powerful Laser source of 30 watts

(Cat No.:33004)

- A combination of three wavelengths: 810 nm, 980 nm & 1064 nm.
- Wavelength mixture: The user can select the % of each wavelength depending on desired clinical effect.
- Automatic selection of suitable wavelength on choosing the required effect.
- Continuous
- Pulsed: Frequency Variable up to 2000Hz
- Duty Cycle variable up to 50%
- High-Intensity Pulse: Sharp pulses for treatment at maximum power and depth with adjustable heat feeling
- Automatic selection of suitable wavelength on choosing the required effect.
- Inbuilt Pathology Library
- Large 10-inch Colour Touch Screen



Contact Us

022-28947733

sales@biotechindia.net

www.biotechstriker.com

